URGENT FIELD SAFETY NOTICE

Potential Cardiac Surgery Infection Risk Reduction Measure: 3T Heater Cooler Design & Operating Instructions update

FSCA Identifier: 9611109-05/05/17-004-C & CP-MUN-2018-009

Affected Devices: All 3T Heater-Cooler Devices

Date: March 20, 2019

Attention: Users of the heater cooler devices (Cardiac Surgery Operating Room Responsible, Cardiac surgeons, Perfusionists, biomedical engineers/technicians, Hygiene Specialists, Risk/ Safety Managers)

Dear Valued Customer:

Reason for this Letter and affected products

The purpose of this letter is to inform you that LivaNova¹ is initiating a design upgrade for 3T heater cooler devices² (3T) in use at your facility.

LivaNova is also releasing an updated version of the Operating Instructions for these devices.

This communication applies to the following Part Numbers:

Part Number	Description
16-02-95	S5 System (Heater-Cooler System 3T)

Description of the addressed issues and corrective measures

1. Design upgrade

In 2017, LivaNova has developed a vacuum aerosol canister and internal sealing design change that is intended to further mitigate (but does not eliminate) the risk of airborne transmission of non-tuberculosis mycobacterium ("NTM") from the 3T device.

The vacuum and sealing solution consists of:

- The modification of the tank lid to improve air tightness, and thus ensure a sufficient level of negative pressure inside the tank when connected to the vacuum source.
- The substitution of the current overflow bottle with a disposable canister connected to a vacuum source. In addition to the collection of water overflowing from the tanks, the canister will collect aerosols. The isolation and collection of the aerosols is achieved by the bacterial filter included in the lid of the vacuum canister to filter the downstream flow prior to the vacuum source.
- The addition of a pressure gauge to measure the vacuum level inside the tank.

¹ LivaNova PLC is a U.K. holding company with a number of wholly-owned subsidiaries, including LivaNova Deutschland GmbH. In this document, we refer to all entities using the brand name LivaNova.

² The 3T is a non sterile heating-cooling machine manufactured by LivaNova Deutschland GmbH. It is used to control patient's body temperature over a cardiopulmonary bypass procedure.

As part of our continuous improvement process, additional design changes aiming to address specific issues, identified through our post-market surveillance system, have also been developed:

- Evaporator coil erosion: the current design of the stirrer pump inside the Patient tank directs water flow to a narrow area of the evaporator cooling coil, thus exacerbating the erosion of the surface of the cooling coil in this location. Ultimately, erosion can result in the creation of a hole in the evaporator, affecting cooling function performance, and in some cases, triggering total shutdown of the device. The newly designed stirrer pump will redirect this water flow away from the cooling coil reducing the potential for erosion on the surface of the cooling coil in this location.
- Change over valve³ failure or performances issues: the current design of the change over valve block for the Cardioplegia tanks allows humid air charged with chemicals used during disinfection to enter into the stepper motor, thus leading to the corrosion of a stepper motor bearing that controls the selection of the Warm or Cold Cardioplegia tank. Change over valve failure can lead to the cardioplegia solution not being cooled or warmed. The design change will prevent humid air from entering the valve block assembly.
- Air vent valves (AVV)⁴ leakage: corrosion residues may hinder the AVV to close properly, resulting in leakages in certain positions of the ball valve. The new AVV design will allow yearly replacement of this part during preventive maintenance operations. This event is not likely to cause any harm to the patient or user.
- Water drain valve⁵ leakage: investigation revealed material abrasion of an internal component in the valve (called the O-ring). A highly abraded O-ring may cause leakage from the valve in the valve open position. The O-ring will be replaced during preventive maintenances operations. This event is not likely to cause any harm to the patient or user.

2. Operating Instructions update

Since 2015, LivaNova has issued different Field Safety Notices to provide additional instructions for the users in relation to potential Nontuberculous Mycobacteria (NTM) contamination Risk.

LivaNova is now releasing a new revision of the Operating Instructions that is applicable to all 3T devices already equipped with the vacuum and sealing design solution. The new Operating Instructions will not directly affect the patient but will consolidate the content of the previously released Field Safety Notices.

What actions should be taken by the Customer/User?

LivaNova will replace all the 3T in use at your facility with newly manufactured devices, incorporating all design changes described in this communication. Replacement devices will be provided to your facility with the new revision of the Operating Instructions. The previous version of the Operating Instructions will not be applicable anymore: it shall remain with replaced 3T or be scrapped.

Your Livanova representative will contact you in the coming weeks to organize with you the replacement of all 3T in use at your facility.

Transmission of this Field Safety Notice

Please ensure that this Field Safety Notice is passed on to all personnel within your organization who need to be aware of it.

In case you have transferred products to a third party, please pass this information on to them and also inform the below mentioned contact person.

³ The change over valve allows the user to switch between cold and warm Cardioplegia tank of the 3T to get the desired temperature range.

⁴ These 3 valves are used to open and close the corresponding water circuit.

⁵ This valve is used to drain the heater cooler device.

Health innovation that matters

Contact reference person

For questions regarding this Field Safety Notice, please contact your LivaNova representative [local contact information to be inserted], or LivaNova Customer Quality at LivaNova.FSCA@livanova.com.

A copy of this Field Safety Notice has been provided to the appropriate Regulatory Agency in your country who is aware of these actions.

Thank you for your cooperation in this matter. LivaNova is committed to providing quality products and we apologize for any inconvenience this may have caused.

Sincerely,

A handwritten signature in blue ink that reads "Galasso".

Barbara Galasso

Customer Quality Manager, Cardiac Surgery

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