


D903 ADVANT

ADULT HOLLOW FIBER OXYGENATOR




SUCTION BLOOD MANAGEMENT

Suction blood is one of the most relevant issue in ECC today.

Suction blood is now recognised as the first culprit in plasma free haemoglobin and platelet activation.

Activated coagulation factors as well as activated mediators are largely present in suction blood.

It has been demonstrated that the separation and the washing treatment of suction blood is of real benefit to the patient. Avant sets a new standard by giving you this choice at no extra cost.



MONITOR PRESSURE

When suction blood is mixed with venous blood and then reinfused into the patient after being oxygenated, the patient's mean arterial pressure drops by 40%. The original value is then re-established over about 10 minutes. 200 mls of reinfused blood is enough to note this phenomenon.

Haemolysis and coagulation system activation is greatly reduced by means of suction blood separation


Exclusive dual reservoir with separate management.

Simple circuit set-up meets any operating theatre configuration.

Reduced tubing lengths.

True management of suctioned blood.

Perfect management of air emboli.



REFERENCES

Perfusion 2001; 16: 519-524

Haemolysis during cardiopulmonary bypass: how to reduce the free haemoglobin by managing the suctioned blood separately

A Pierangeli et Al - St. Orsola University Hospital, Bologna

The Journal of Thoracic and Cardiovascular Surgery - May 2002

Tissue factor as the main activator of the coagulation system during cardiopulmonary bypass

De Somer et Al - University Hospital Gent

When you begin an ECC you never know how much suction blood you will get back. Avant gives you the possibility of deciding at any time whether it is wiser to reject the blood or wash it, at no extra work.




Field suction

Venous return

Left vent


Arterial filter purge

Haemofilter line



The flow between the cardiotomy and venous reservoirs is controlled by a blue lever.

The connection is fully open when the key is up and closed when down.




"Aspirated blood contaminated by tissue contact is the most important activator of the coagulation system and the principal cause of hemolysis during cardiopulmonary bypass."
(J Thorac Cardioasc Surg 2002; 123:951-8)

SAFETY AND FLEXIBILITY

SAFE.


The safety of the patient is fundamental during bypass.



The contribution of the suckers to bubble activity is evident. Avant has a separate cardiomy reservoir which helps remove bubbles before the blood enters the venous reservoir.

In an open system the risk of emboli to the patient due to a low level within the reservoir is often reduced by a low level detector. When, however, this system is not available or fails, Sorin Group patented safety valve comes to your help.

The clinical experience of this valve over 1.500.000 cases during the last decade, is proof of its value. Past and present users confirm the absence of valve collapse and the benefits to patient safety provided by the safety valve.



FLEXIBLE.


Flexibility of the system, together with the Avant holder, allow you to move the oxygenator in any position and to reduce tubing length to easily adapt to any operating theatre configuration.



The 180° rotatable venous return connector gives access and visibility to the venous line whichever side you want it.



The holder has integrated Hansen water connectors which secure Avant when the gold indicates "closed".




The 250° independent rotation of the oxy module and hard shell venous reservoir matches any pump configuration. Left or right handed devices are now a thing of the past.



The 360° rotatable turrets help shorten suckers and therefore simplify the circuit.



When Avant is set on the holder it is possible to independently rotate the oxy module and venous reservoir as well as the whole unit.




The holder has been designed to make life easier.

COATING

Once the suctions have been managed, the patient can truly benefit from a non-thrombogenic surface in ECC. Ph.I.S.I.O. is a second generation, truly biomimetic coating which aims at improving the haemocompatibility of surfaces in contact with blood

Thanks to the presence of phosphorylcholine polar groups in the Ph.I.S.I.O. polymer there is a phospholipid-like structure that creates a biomembrane on the surfaces when they come in contact with the blood cell membrane phospholipids.


The polymer's physical structure consists of hydrophilic and hydrophobic parts and contains balanced presence of positive and negative charges. In this way it performs as an optimal interface between polymeric surfaces and blood.



Reduced platelets adhesion.


Only one kind of surface in contact with blood.

Phisio coating and Avant suction separation held in reducing post-op. bleeding. (Perfusion 2002; 17:39-44)



REFERENCES





- Perfusion 2002; 17: 39-44
Phosphorylcholine coating offers natural platelet preservation during cardiopulmonary bypass
De Somer et Al - Gent University Hospital
- The Journal of Thoracic and Cardiovascular Surgery - November 2003
Physiological coagulation can be maintained in extracorporeal circulation by means of shed blood separation and coating
Albes et Al - Gena University Hospital
- European Journal of Cardio-thoracic surgery 18 (2000) 602-606
Phosphorylcholine coating of extracorporeal circuits provides natural protection against blood activation by the material surface
De Somer et Al - Gent University Hospital




Phisio has proven to be extremely effective in reducing platelet activation and platelet foreign surface adhesion.

"Combination of coating and avoidance of shed blood recirculation maintained physiological coagulation levels and markedly reduced red blood cell trauma in extracorporeal circulation procedures. These combined modalities may therefore offer an alternative for off-pump procedures in patients with contraindications for conventional extracorporeal circulation"
(J Thorac Cardiovasc Surg 2003; 126:1504-12)

PERFORMANCE



TECHNICAL FEATURES



Max
Blood Flow
7.5 l/min

Static
Priming Volume
(recovered)
250 ml

Membrane
Surface Area
1.7 sq.m.

H.E. Performance
Factor
h 0.68 (4 l/min)

Cadiotomy Filter
Pore Size
30 um

オーダーガイド

製品番号 型番 製品名

05332J D903 D903アヴァント フィジオ*(人工肺・リザーバー 一体型、セーフティバルブ付)

05331J D903 D903アヴァント フィジオ*(人工肺モジュール)

*フィジオ:リン脂質ポリマーコーティング『フィジオ』付

05365J D970 D903アヴァント(ハードシェル静脈リザーバー単体)

ホルダー:

05046 D632 成人肺一体型ホルダー(熱交換器カブラー付)

05047 - 成人肺一体型ホルダー(熱交換器カブラー無)

050128 - 成人肺単体用ホルダー


05042 D628 D903アヴァント ハードシェル静脈リザーバー単体用ホルダー

付属品:

05490J D890 サンプルングマニホールド

05069 D712 サンプルングマニホールド用ホルダー

販 売 名	医療機器承認番号
D903アヴァント	20900BZY00560000
D903アヴァントPhisio	21700BZG00006000
選 任 製 造 販 売 業 者	ソーリン・グループ株式会社 〒100-6110東京都千代田区永田町2-11-1 Tel:03-3595-7630(代) Fax:03-3595-7631
外国特例承認取得者	ソーリン・グループ・イタリア社:イタリア



The Sorin Group Italia Quality
System complies with:
EN ISO 13485:2003/AC:2007

MANUFACTURED BY:
SORIN GROUP ITALIA S.R.L.
VIA STATALE 12 NORD, 86
41037 MIRANDOLA MODENA ITALY
TEL +39 0535 29811
FAX +39 0535 24312
INFO.CP@SORIN.COM

WWW.SORIN.COM