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
**Monitor**  
**Barom.**  
**Trends**  
**Set-Up**

**Alarms**  
**Time/Date**  
**RS-232**

Calibrating venous probe pls. wait 042sec.

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Arterial cable connected. pO<sub>2</sub> sensor not connected. Check connection and wait a few seconds. To continue with venous monitoring only press ENTER.

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Arterial cable connected.  
pO<sub>2</sub> sensor connected.  
pO<sub>2</sub> calibration underway.

Venous cable connected. Calibration on stand-by. Clean reference cell and correctly park probe onto it. Check and wait 15 seconds. To continue without venous monitoring press ENTER.

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| Guidance and manufacturer's declaration – electromagnetic emission – for all EQUIPMENT AND SYSTEMS (see 5.2.2.1 c))                                                                       |                |                                                                                                                                                                                                           |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Guidance and manufacturer's declaration – electromagnetic emission                                                                                                                        |                |                                                                                                                                                                                                           |
| The DATA MASTER is intended for use in the electromagnetic environment specified below. The customer or the user of the DATA MASTER should assure that it is used in such an environment. |                |                                                                                                                                                                                                           |
| Emissions test                                                                                                                                                                            | Compliance     | Electromagnetic environment – guidance                                                                                                                                                                    |
| RF emissions<br>CISPR 11                                                                                                                                                                  | <b>Group 1</b> | The DATA MASTER uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.                      |
| RF emissions<br>CISPR 11                                                                                                                                                                  | Class A        | The DATA MASTER is suitable for use in all establishments other than domestic and those directly connected to the public-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions<br>IEC 61000-3-2                                                                                                                                                       | Class A        |                                                                                                                                                                                                           |
| Voltage fluctuations / flicker emissions<br>IEC 61000-3-3                                                                                                                                 | Complies       |                                                                                                                                                                                                           |

| <b>Guidance and manufacturer's declaration – electromagnetic immunity – for all me equipment and systems (see 5.2.2.1 f))</b>                                                                |                                                                                                                                                                                                             |                                                                                                                                                                                                             |                                                                                                                                                                                                                                                                                        |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Guidance and manufacturer's declaration – electromagnetic immunity</b>                                                                                                                    |                                                                                                                                                                                                             |                                                                                                                                                                                                             |                                                                                                                                                                                                                                                                                        |
| The DATA MASTER is intended for use in the electromagnetic environment specified below.<br>The customer or the user of the DATA MASTER should assure that it is used in such an environment. |                                                                                                                                                                                                             |                                                                                                                                                                                                             |                                                                                                                                                                                                                                                                                        |
| <b>IMMUNITY TEST</b>                                                                                                                                                                         | <b>IEC 60601 test level</b>                                                                                                                                                                                 | <b>Compliance level</b>                                                                                                                                                                                     | <b>Electromagnetic environment - guidance</b>                                                                                                                                                                                                                                          |
| Electrostatic discharge (ESD)<br><br>IEC 61000-4-2                                                                                                                                           | ± 6 kV contact<br><br>± 8 kV air                                                                                                                                                                            | ± 6 kV contact<br><br>± 8 kV air                                                                                                                                                                            | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.                                                                                                                                         |
| Electrostatic transient / burst<br><br>IEC 61000-4-4                                                                                                                                         | ± 2 kV for power supply lines<br>± 1 kV for input/output lines                                                                                                                                              | ± 2 kV for power supply lines<br>± 1 kV for input/output lines                                                                                                                                              | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                                    |
| Surge<br><br>IEC 61000-4-5                                                                                                                                                                   | ± 1 kV differential mode<br>± 2 kV common mode                                                                                                                                                              | ± 1 kV differential mode<br>± 2 kV common mode                                                                                                                                                              | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                                    |
| Voltage dips, short interruptions and voltage variations on power supply input lines<br><br>IEC 61000-4-11                                                                                   | < 5 % $U_T$<br>(>95 % dip in $U_T$ ) for 0,5 cycle<br>40 % $U_T$<br>(60 % dip in $U_T$ ) for 5 cycles<br>70 % $U_T$<br>(30 % dip in $U_T$ ) for 25 cycles<br>< 5 % $U_T$<br>(>95 % dip in $U_T$ ) for 5 sec | < 5 % $U_T$<br>(>95 % dip in $U_T$ ) for 0,5 cycle<br>40 % $U_T$<br>(60 % dip in $U_T$ ) for 5 cycles<br>70 % $U_T$<br>(30 % dip in $U_T$ ) for 25 cycles<br>< 5 % $U_T$<br>(>95 % dip in $U_T$ ) for 5 sec | Mains power quality should be that of a typical commercial or hospital environment. If the user of the DATA MASTER requires continued operation during power mains interruptions, it is recommended that the DATA MASTER be powered from an uninterruptible power supply or a battery. |
| Power frequency (50/60 Hz) magnetic field<br><br>IEC 61000-4-8                                                                                                                               | 3 A/m                                                                                                                                                                                                       | 3 A/m                                                                                                                                                                                                       | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.                                                                                                                                              |
| NOTE $U_T$ is the a. c. mains voltage prior to application of the test level.                                                                                                                |                                                                                                                                                                                                             |                                                                                                                                                                                                             |                                                                                                                                                                                                                                                                                        |

| Guidance and MANUFACTURER’S declaration – ELECTROMAGNETIC IMMUNITY – for ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING (see 5.2.2.2)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |                             |                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------|------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Guidance and manufacturer’s declaration – electromagnetic immunity                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |                             |                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| The DATA MASTER is intended for use in the electromagnetic environment specified below. The customer or the user of the DATA MASTER should assure that it is used in such an environment.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                             |                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| IMMUNITY TEST                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       | IEC 60601 test level        | Compliance level | Electromagnetic environment - guidance                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Conducted RF<br>IEC 61000-4-6                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       | 3 Vrms<br>150 kHz to 80 MHz | 3 Vrms           | Portable and mobile RF communications equipment should be used no closer to any part of the DATA MASTER, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.<br><br><b>RECOMMENDED SEPARATION DISTANCE</b><br>$d = 1,2\sqrt{P}$                                                                                                                                                                                                                                                                                                                                                              |
| Conducted RF<br>IEC 61000-4-3                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       | 3 V/m<br>80 MHz to 2.5 GHz  | 3 V/m            | $d = 1,2\sqrt{P}$ 80 MHz to 800 MHz<br><br>$d = 2,3\sqrt{P}$ 800 MHz to 2,5 GHz<br><br>where $p$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and $d$ is the recommended separation distance in metres (m). <sup>b</sup><br><br>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup><br><br>Interference may occur in the vicinity of equipment marked with the following symbol:<br><br> |
| NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |                             |                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| NOTE 2 These guidelines may not apply in all situations. Electromagnetic is affected by absorption and reflection from structures, objects and people.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |                             |                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the DATA MASTER is used exceeds the applicable RF compliance level above, the DATA MASTER should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the DATA MASTER. |                             |                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |                             |                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |

60601 – 1 – 2/FDIS © IEC(E)

**and the ME EQUIPMENT or ME SYSTEM – for ME EQUIPMENT and ME SYSTEMS that are not LIFE SUPPORTING (see 5.2.2.2)**

**Recommended separation distances between portable and mobile RF communications equipment and the DATA MASTER**

The DATA MASTER is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the DATA MASTER can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the DATA MASTER as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum output of transmitter<br><br>W | Separation distance according to frequency of transmitter<br>M |                                              |                                        |
|----------------------------------------------|----------------------------------------------------------------|----------------------------------------------|----------------------------------------|
|                                              | 150 kHz to 80 MHz<br>$d = 1,2\sqrt{P}$                         | Rated maximum output of transmitter<br><br>W | 150 kHz to 80 MHz<br>$d = 1,2\sqrt{P}$ |
| 0,01                                         | 0,12                                                           | 0,01                                         | 0,12                                   |
| 0,1                                          | 0,38                                                           | 0,1                                          | 0,38                                   |
| 1                                            | 1,2                                                            | 1                                            | 1,2                                    |
| 10                                           | 3,8                                                            | 10                                           | 3,8                                    |
| 100                                          | 12                                                             | 100                                          | 12                                     |

For transmitters rated at a maximum output power not listed above the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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